

Anti-FGF4 Antibody Biotin Conjugated

Catalog Number: A01675-1-Biotin

About FGF4

Fibroblast growth factor 4(FGF4), also known as Heparin Secretary Transforming(HSTF1). HST1, for which the designation HSTF1 was proposed for human gene nomenclature, is a heparin-binding growth factor with significant homology to human fibroblast growth factors and the mouse Int-2 protein. By in situ hybridization, Adelaide et al.(1988) mapped the HST gene to chromosome 11q13. The HST1 protein is a heparin-binding growth factor with significant homology with human fibroblast growth factors and the mouse Int-2 protein.

Overview

Product Name	Anti-FGF4 Antibody Biotin Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P08620

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human FGF4.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

1 Publications Citing This Product

1. PubMed ID: 10.1016/j.bbrc.2016.01.174, The oncoprotein HBXIP up-regulates FGF4 through activating transcriptional factor Sp1 to promote the migration of breast cancer cells

Visit bosterbio.com/anti-fgf4-picoband-trade-antibody-a01675-1-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FGF4 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.