

Anti-Human CCKBR DyLight® 550 conjugated Antibody

Catalog Number: A01677-DyI550

About CCKBR

The cholecystokinin B receptor, also known as CCKBR or CCK2, is a protein that in humans is encoded by the CCKBR gene. This gene encodes a G-protein coupled receptor for gastrin and cholecystokinin (CCK), regulatory peptides of the brain and gastrointestinal tract. This protein is a type B gastrin receptor, which has a high affinity for both sulfated and nonsulfated CCK analogs and is found principally in the central nervous system and the gastrointestinal tract. Alternative splicing results in multiple transcript variants. A misspliced transcript variant including an intron has been observed in cells from colorectal and pancreatic tumors.

Overview

Product Name	Anti-Human CCKBR DyLight® 550 conjugated Antibody
Reactive Species	Human
Description	Boster Bio Anti-Human CCKBR DyLight® 550 conjugated Antibody catalog # A01677-DyI550. Tested in Flow Cytometry applications. This antibody reacts with Human.
Conjugate	DyLight®550
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P32239

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human CCKBR, which shares 80% and 90% amino acid (aa) sequence identity with mouse and rat CCKBR, respectively.
Predicted Reactive Species	Human
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.

Suggested Dilutions

Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.

If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.

Some PubMed article(s) citing the expression level of this target are as follows:

Boster Bio's internal QC testing used:

Flow Cytometry, 1-3ug/1x10⁶ cells

Anti-Human CCKBR DyLight® 550 conjugated Antibody (A01677-DyI550) Images



Boster Kit Box

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Human CCKBR DyLight® 550 conjugated Antibody