

Anti-RFC/SLC19A1 Antibody Picoband® HRP Conjugated

Catalog Number: A01679-1-HRP

About SLC19A1

Solute carrier family 19 (folate transporter), member 1, also known as SLC19A1 or RFC1, is a protein which in humans is encoded by the SLC19A1 gene. Transport of folate compounds into mammalian cells can occur via receptor-mediated or carrier-mediated mechanisms. A functional coordination between these 2 mechanisms has been proposed to be the method of folate uptake in certain cell types. Methotrexate (MTX) is an antifolate chemotherapeutic agent that is actively transported by the carrier-mediated uptake system. And RFC1 plays a role in maintaining intracellular concentrations of folate.

Overview

Product Name	Anti-RFC/SLC19A1 Antibody Picoband® HRP Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P41440

Technical Details

Immunogen	E.coli-derived human RFC/SLC19A1 recombinant protein (Position: V198-N231).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RFC/SLC19A1 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.