

## Anti-DNA Ligase IV LIG4 Antibody

Catalog Number: A01685-2


### Overview

Product Name	Anti-DNA Ligase IV LIG4 Antibody
Reactive Species	Human
Description	Boster Bio Anti-DNA Ligase IV LIG4 Antibody catalog # A01685-2. Tested in WB,IHC,ICC applications. This antibody reacts with Human.
Application	IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P49917


### Technical Details

Immunogen	Recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 IHC: 1:50-1:200 ICC: 1:10-1:50

## Anti-DNA Ligase IV LIG4 Antibody (A01685-2) Images



Western blot analysis of DNA Ligase IV on HeLa cell using anti-DNA Ligase IV antibody at 1/1,000 dilution. Positive control: Lane 1: HepG2 Lane 2: HeLa Lane 3: Jurkat



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue using anti-DNA Ligase IV antibody. Counter stained with hematoxylin.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DNA Ligase IV LIG4 Antibody

For Research Use Only. Not for use in diagnostic procedures.