

Anti-CD32-B FCGR2B Antibody

Catalog Number: A01690


Overview

| | |
|----------------------|--|
| Product Name | Anti-CD32-B FCGR2B Antibody |
| Reactive Species | Human |
| Description | Boster Bio Anti-CD32-B FCGR2B Antibody catalog # A01690. Tested in WB applications. This antibody reacts with Human. |
| Application | WB |
| Clonality | Polyclonal |
| Formulation | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | P31994 |

Technical Details

| | |
|---------------------|---|
| Immunogen | Recombinant fusion protein of human FCGR2B(NP_003992.3). |
| Cross Reactivity | No cross reactivity with other proteins. |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | WB: 1:500-1:2000 |

Anti-CD32-B FCGR2B Antibody (A01690) Images



Western blot analysis of extracts of human placenta, using FCGR2B antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD32-B FCGR2B Antibody

For Research Use Only. Not for use in diagnostic procedures.