

Anti-CD32-B (phospho-Y292) FCGR2B Antibody

Catalog Number: A01690Y292


Overview

Product Name	Anti-CD32-B (phospho-Y292) FCGR2B Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-CD32-B (phospho-Y292) FCGR2B Antibody catalog # A01690Y292. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P31994

Technical Details

Immunogen	Synthetic phosphopeptide derived from human CD32-B around the phosphorylation site of Tyrosine 292.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-CD32-B (phospho-Y292) FCGR2B Antibody (A01690Y292) Images



Western blot (WB) analysis of p-CD32-B (Y292) polyclonal antibody at 1:500 dilution Lane1:MCF-7 cell lysate treated with pervanadate (25ng/ML,72h) Lane2:Mouse liver tissue lysate Lane3:PC12 cell lysate treated with pervanadate (25ng/ML,72h)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CD32-B (phospho-Y292) FCGR2B Antibody

For Research Use Only. Not for use in diagnostic procedures.