

## Anti-40S ribosomal protein SA RPSA Antibody

Catalog Number: A01691-2


### Overview

|                      |   |
|----------------------|---|
| Product Name         | Anti-40S ribosomal protein SA RPSA Antibody   |
| Reactive Species     | Human, Mouse, Rat   |
| Description          | Boster Bio Anti-40S ribosomal protein SA RPSA Antibody catalog # A01691-2. Tested in WB,ICC/IF,IHC,IP,Flow Cytometry applications. This antibody reacts with Human,Mouse,Rat. |
| Application          | Flow Cytometry, IP, IF, IHC, ICC, WB  |
| Clonality            | Polyclonal  |
| Formulation          | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2  |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.                                    |
| Host                 | Rabbit  |
| Uniprot ID           | P08865  |


### Technical Details

|                     |  |
|---------------------|--|
| Immunogen           | recombinant protein  |
| Cross Reactivity    | No cross reactivity with other proteins.                                     |
| Isotype             | IgG  |
| Form                | Liquid   |
| Concentration       | 1 mg/ml  |
| Purification        | ProA affinity purified   |
| Suggested Dilutions | WB: 1:1,000-1:2,000<br>ICC: 1:100-1:500<br>IHC: 1:50-1:200<br>FC: 1:50-1:100 |

## Anti-40S ribosomal protein SA RPSA Antibody (A01691-2) Images



Western blot analysis of RPSA on different lysates using anti-RPSA antibody at 1/1,000 dilution. Positive control: Lane 1: K562 Lane 2: HepG2 Lane 3: A431



ICC staining RPSA in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-40S ribosomal protein SA RPSA Antibody

For Research Use Only. Not for use in diagnostic procedures.