

Anti-DDR2/Tyro10 Antibody

Catalog Number: A01698


Overview

Product Name	Anti-DDR2/Tyro10 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-DDR2/Tyro10 Antibody catalog # A01698. Tested in IF, IHC, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, IHC
Clonality	Polyclonal
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q16832


Technical Details

Immunogen	Synthetic peptide from human protein at AA range: 31-80
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	IHC 1:50-200 IF 1:500 ELISA 1:10000-20000

Anti-DDR2/Tyro10 Antibody (A01698) Images



Immunohistochemical analysis of paraffin-embedded human-lung-cancer, antibody was diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded human-lung-cancer, antibody was diluted at 1:200.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DDR2/Tyro10 Antibody

For Research Use Only. Not for use in diagnostic procedures.