

Anti-Troponin I3 (phospho-S22/S23) TNNI3 Antibody

Catalog Number: A01720S22S23


Overview

| | |
|----------------------|--|
| Product Name | Anti-Troponin I3 (phospho-S22/S23) TNNI3 Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-Troponin I3 (phospho-S22/S23) TNNI3 Antibody catalog # A01720S22S23. Tested in WB applications. This antibody reacts with Human, Mouse, Rat. |
| Application | WB |
| Clonality | Polyclonal |
| Formulation | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | P19429 |

Technical Details

| | |
|---------------------|---|
| Immunogen | Synthetic phosphopeptide derived from human Troponin I3 around the phosphorylation site of Ser22/23. |
| Cross Reactivity | No cross reactivity with other proteins. |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | WB: 1:500-1:1000 |

Anti-Troponin I3 (phospho-S22/S23) TNNI3 Antibody (A01720S22S23) Images



Western blot (WB) analysis of p-Troponin I3 (S22/S23) pAb at 1:500 dilution Lane1:The Heart tissue lysate of Mouse(10ug) Lane2:The Heart tissue lysate of Rat(10ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Troponin I3 (phospho-S22/S23) TNNI3 Antibody

For Research Use Only. Not for use in diagnostic procedures.