

Anti-PAK4/5/6 (phospho-S474) Antibody

Catalog Number: A01723S474


Overview

Product Name	Anti-PAK4/5/6 (phospho-S474) Antibody
Reactive Species	Human
Description	Boster Bio Anti-PAK4/5/6 (phospho-S474) Antibody catalog # A01723S474. Tested in WB,IHC,ICC applications. This antibody reacts with Human.
Application	IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O96013/Q9NQ5/Q9P286

Technical Details

Immunogen	Synthetic phosphopeptide derived from human PAK4 around the phosphorylation site of Serine 474.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 ICC: 1:50-1:200

Anti-PAK4/5/6 (phospho-S474) Antibody (A01723S474) Images



Western blot (WB) analysis of p-PAK4/5/6 (S474) pAb at 1:500 dilution Lane1:PC3 whole cell lysate(40ug) Lane2:HCT116 whole cell lysate(40ug) Lane3:MCF-7 whole cell lysate(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PAK4/5/6 (phospho-S474) Antibody

For Research Use Only. Not for use in diagnostic procedures.