

Anti-MYPT1 (phospho-T696) PPP1R12A Antibody

Catalog Number: A01743T696


Overview

Product Name	Anti-MYPT1 (phospho-T696) PPP1R12A Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-MYPT1 (phospho-T696) PPP1R12A Antibody catalog # A01743T696. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O14974

Technical Details

Immunogen	Synthetic peptide derived from human MYPT1 around the phosphorylation site of T696.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-MYPT1 (phospho-T696) PPP1R12A Antibody (A01743T696) Images



Western blot (WB) analysis of p-MYPT1 (T696) polyclonal antibody at 1:500 dilution Lane1:HEK293T whole cell lysate(40ug) Lane2:HEK293T treated with PMA (60ng/ml, 10min) whole cell lysate (40ug) Lane3:HEK293T treated with UV (4h) whole cell lysate (40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MYPT1 (phospho-T696) PPP1R12A Antibody

For Research Use Only. Not for use in diagnostic procedures.