

## Anti-MAP1B Antibody

Catalog Number: A01760

### Overview

Product Name	Anti-MAP1B Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-MAP1B Antibody catalog # A01760. Tested in IHC, IF applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, IHC
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P46821

### Technical Details

Immunogen	Synthesized peptide derived from human protein at AA range: 480-560
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	MAP1B antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	IHC 1:50-300 IF 1:50-200

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MAP1B Antibody

For Research Use Only. Not for use in diagnostic procedures.