

Anti-PKC Zeta (phospho-T410) PRKCZ Antibody

Catalog Number: A01796T410


Overview

Product Name	Anti-PKC Zeta (phospho-T410) PRKCZ Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PKC Zeta (phospho-T410) PRKCZ Antibody catalog # A01796T410. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q05513


Technical Details

Immunogen	Synthetic phosphopeptide derived from human PKC zeta around the phosphorylation site of Threonine 410.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IF: 1:50-1:200

Anti-PKC Zeta (phospho-T410) PRKCZ Antibody (A01796T410) Images



Western blot (WB) analysis of p-PKC zeta (T410)pAb at 1:500 dilution Lane1:MCF-7 whole cell lysate Lane2:MCF-7 treated with insulin(100nM/ml,15 minutes) whole cell lysate Lane3:The Kidney tissue lysate of Rat Lane4:The Lung tissue lysate of Mouse



Immunohistochemistry (IHC) analyzes of p-PKC zeta (T410) pAb in paraffin-embedded human breast carcinoma tissue at 1:100.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PKC Zeta (phospho-T410) PRKCZ Antibody

For Research Use Only. Not for use in diagnostic procedures.