

## Anti-PKC Zeta (phospho-T560) PRKCZ Antibody

Catalog Number: A01796T560


### Overview

Product Name	Anti-PKC Zeta (phospho-T560) PRKCZ Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PKC Zeta (phospho-T560) PRKCZ Antibody catalog # A01796T560. Tested in WB,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q05513


### Technical Details

Immunogen	Synthetic phosphopeptide derived from human PKC zeta around the phosphorylation site of Threonine 560.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IP: 1:10-1:100


## Anti-PKC Zeta (phospho-T560) PRKCZ Antibody (A01796T560) Images



Western blot (WB) analysis of p-PKC zeta (T560) pAb at 1:500 dilution Lane1:The Brain tissue lysate of Mouse(40ug) Lane2:The Brain tissue lysate of Rat(40ug) Lane3:HEK293T whole cell lysate(20ug) Lane4:PC3 whole cell lysate(20ug)



Immunohistochemistry (IHC) analyzes of p-PKC zeta (T560) pAb in paraffin-embedded human brain tissue.



Immunoprecipitation of Brain tissue lysate of Rat using p-PKC zeta (T560) polyclonal antibody (Sepharose Bead Conjugate) #BD0048(lane 2 and lane 3) .Lane 1 is 30% input.The western blot was probed using p-PKC zeta (T560) #BS4347.“↑”supernatant;“↓”(deposition)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PKC Zeta (phospho-T560) PRKCZ Antibody

For Research Use Only. Not for use in diagnostic procedures.