

Anti-NMDA Zeta 1 (phospho-S896) GRIN1 Antibody

Catalog Number: A01808S896


Overview

Product Name	Anti-NMDA Zeta 1 (phospho-S896) GRIN1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-NMDA Zeta 1 (phospho-S896) GRIN1 Antibody catalog # A01808S896. Tested in IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q05586

Technical Details

Immunogen	Synthetic phosphopeptide derived from human NMDAzeta1 around the phosphorylation site of Serine 896.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200

Anti-NMDA Zeta 1 (phospho-S896) GRIN1 Antibody (A01808S896) Images



Immunohistochemistry (IHC) analyzes of p-NMDAZeta1 (S896) pAb in paraffin-embedded human breast carcinoma tissue at 1:100.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NMDA Zeta 1 (phospho-S896) GRIN1 Antibody

For Research Use Only. Not for use in diagnostic procedures.