

## Anti-ERM/Etv5 Antibody Picoband® Biotin Conjugated

Catalog Number: A01809-2-Biotin

### About ETV5

Ets variant 5 (ETV5) (also named ERM transcription factor) is a transcription factor that in humans is encoded by the ETV5 gene. It is generated in Sertoli cells, which are found in the testes and play a crucial role in spermatogenesis. Its ortholog has been linked to both obesity and bipolar disorder.

### Overview

Product Name	Anti-ERM/Etv5 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P41161

### Technical Details

Immunogen	E.coli-derived human ERM/Etv5 recombinant protein (Position: D72-Q246).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-ERM/Etv5 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.