

## Anti-JunB (phospho-S79) Antibody

Catalog Number: A01825S79

### Overview


Product Name	Anti-JunB (phospho-S79) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-JunB (phospho-S79) Antibody catalog # A01825S79. Tested in IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P17275

### Technical Details

Immunogen	Synthetic phosphopeptide derived from human JunB around the phosphorylation site of Serine 79.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200 IP: 1:50-1:200

## Anti-JunB (phospho-S79) Antibody (A01825S79) Images

---



Immunohistochemistry (IHC) analysis of p-JunB (S79) pAb in paraffin-embedded human breast carcinoma tissue.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

---

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-JunB (phospho-S79) Antibody

For Research Use Only. Not for use in diagnostic procedures.