

Anti-ETAR EDNRA Antibody

Catalog Number: A01828


Overview

| | |
|----------------------|--|
| Product Name | Anti-ETAR EDNRA Antibody |
| Reactive Species | Human, Monkey, Mouse |
| Description | Boster Bio Anti-ETAR EDNRA Antibody catalog # A01828. Tested in WB applications. This antibody reacts with Human, Monkey, Mouse. |
| Application | WB |
| Clonality | |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | P25101 |

Technical Details

| | |
|---------------------|--|
| Immunogen | Synthesized peptide derived from the C-terminal region of human ETAR. at AA range: 350-430 |
| Isotype | IgG |
| Form | Liquid |
| Purification | Immunogen affinity purified |
| Suggested Dilutions | <p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used WB, 1:500-1:2000</p> <p>For protocols please visit https://www.bosterbio.com/protocol-and-troubleshooting/</p> |

Anti-ETAR EDNRA Antibody (A01828) Images



Western Blot (WB) analysis of specific cells using ETAR Polyclonal antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ETAR EDNRA Antibody