

## Anti-Endothelin A Receptor/ET-A/EDNRA Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A01828-4-Fluoro594

### About EDNRA

This gene encodes the receptor for endothelin-1, a peptide that plays a role in potent and long-lasting vasoconstriction. This receptor associates with guanine-nucleotide-binding (G) proteins, and this coupling activates a phosphatidylinositol-calcium second messenger system. Polymorphisms in this gene have been linked to migraine headache resistance. Alternative splicing results in multiple transcript variants.

### Overview

Product Name	Anti-Endothelin A Receptor/ET-A/EDNRA Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P25101

### Technical Details

Immunogen	E.coli-derived human Endothelin A Receptor/ET-A/EDNRA recombinant protein (Position: N32-Q381).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

## 1 Publications Citing This Product

1. PubMed ID: 20167215, Peng Hj, Dai Dz, Ji H, Dai Y. Eur J Pharmacol. 2010 May 25;634(1-3):101-6. Doi: 10.1016/J.Ejphar.2010.02.001. Epub 2010 Feb 16. The Separate Roles Of Endothelin Receptors Participate In Remodeling Of Matrix Metalloproteinase And Connexin 43 Of Car...

Visit [bosterbio.com/anti-endothelin-a-receptor-et-a-ednra-picoband-trade-antibody-a01828-4-boster.html](http://bosterbio.com/anti-endothelin-a-receptor-et-a-ednra-picoband-trade-antibody-a01828-4-boster.html) to see all 1 publications.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Endothelin A Receptor/ET-A/EDNRA Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.