

Anti-Protoporphyrinogen oxidase PPOX Antibody

Catalog Number: A01856


Overview

Product Name	Anti-Protoporphyrinogen oxidase PPOX Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Protoporphyrinogen oxidase PPOX Antibody catalog # A01856. Tested in WB,IHC applications. This antibody reacts with Human,Mouse.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P50336


Technical Details

Immunogen	Recombinant fusion protein of human PPOX(NP_001116236.1).
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:50-1:200

Anti-Protoporphyrinogen oxidase PPOX Antibody (A01856) Images



Western blot analysis of extracts of various cell lines, using PPOX Antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.



Immunohistochemistry of paraffin-embedded human colon carcinoma using PPOX antibody at dilution of 1:100. Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Protoporphyrinogen oxidase PPOX Antibody

For Research Use Only. Not for use in diagnostic procedures.