

## Anti-TMEM173 Antibody

Catalog Number: A01871-2


### Overview

Product Name	Anti-TMEM173 Antibody
Reactive Species	Human
Description	Boster Bio Anti-TMEM173 Antibody catalog # A01871-2. Tested in WB,Flow Cytometry,ICC/IF applications. This antibody reacts with Human.
Application	Flow Cytometry, IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q86WV6


### Technical Details

Immunogen	Recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 FC: 1:50-1:100 ICC/IF: 1:10-1:50

## Anti-TMEM173 Antibody (A01871-2) Images



Western blot analysis of TMEM173 on U937 (1) and 293 (2) cell lysate using anti-TMEM173 antibody at 1/1,000 dilution.



Flow cytometric analysis of THP-1 cells with TMEM173 antibody at 1/100 dilution (red) compared with an unlabelled control (cells without incubation with primary antibody; black).

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TMEM173 Antibody

For Research Use Only. Not for use in diagnostic procedures.