

Anti-Xanthine Oxidase XDH Antibody

Catalog Number: A01884-1

Overview

Product Name	Anti-Xanthine Oxidase XDH Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Xanthine Oxidase XDH Antibody catalog # A01884-1. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P47989

Technical Details

Immunogen	Recombinant protein within human Xanthine Oxidase aa 100-300.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 IHC: 1:50-1:200

Anti-Xanthine Oxidase XDH Antibody (A01884-1) Images




Figure 2. Immunohistochemistry validation of XDH using Anti-Xanthine Oxidase XDH Antibody (A01884-1).

Immunohistochemical analysis of paraffin-embedded rat liver tissue using anti-Xanthine Oxidase antibody. Counter stained with hematoxylin.
For more protocol information of IHC




Figure 1. Western blotting validation for Anti-Xanthine Oxidase XDH Antibody A01884-1

Western blot analysis of Xanthine Oxidase on human kidney tissue lysate using anti-Xanthine Oxidase antibody at 1/1

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Xanthine Oxidase XDH Antibody