

## Anti-PKC gamma/PRKCG Antibody Picoband® Biotin Conjugated

Catalog Number: A01890-Biotin

### About PRKCG

The gamma isotype of protein kinase C (PKC gamma) is a member of the classical PKC (cPKC) subfamily which is activated by Ca (2+) and diacylglycerol in the presence of phosphatidylserine. Physiologically, PKC gamma is activated by a mechanism coupled with receptor-mediated breakdown of inositol phospholipid as other cPKC isotypes such as PKC alpha and PKC beta. PKC gamma is expressed solely in the brain and spinal cord and its localization is restricted to neurons, while PKC alpha and PKC beta are expressed in many tissues in addition to the brain. Within the brain, PKC gamma is the most abundant in the cerebellum, hippocampus and cerebral cortex, where the existence of neuronal plasticity has been demonstrated. PKC gamma gene is mutated in spinocerebellar ataxia type 14 (SCA14). Verbeek et al. (2005) point out the specific alterations in mutant PKC gamma function that could lead to the selective neuronal degeneration of SCA14.

### Overview

Product Name	Anti-PKC gamma/PRKCG Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P05129

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human PKC gamma/PRKCG, identical to the related mouse and rat sequences.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

### 3 Publications Citing This Product

1. PubMed ID: 10.2147/DDDT.S93570, 5-(Bis(3-(2-hydroxyethyl)-1H-indol-2-yl)methyl)-2-hydroxybenzoic acid (BHIMHA): showing a strategy of designing drug to block lung metastasis of tumors
2. PubMed ID: 26937173, 5-(Bis(3-(2-hydroxyethyl)-1H-indol-2-yl)methyl)-2-hydroxybenzoic acid (BHIMHA): showing a strategy of designing drug to block lung metastasis of tumors
3. PubMed ID: 25160038, Extracts from Salvia-Nelumbinis naturalis alleviate hepatosteatosi via improving hepatic insulin sensitivity

Visit [bosterbio.com/anti-pkc-gamma-prkcg-picoband-trade-antibody-a01890-boster.html](https://bosterbio.com/anti-pkc-gamma-prkcg-picoband-trade-antibody-a01890-boster.html) to see all 3 publications.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PKC gamma/PRKCG Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.