

Anti-RAB5A/B/C Antibody Picoband® Biotin Conjugated

Catalog Number: A01891-Biotin

About RAB5A/B/C

Members of the Rab protein family are small GTPases of the Ras superfamily that are thought to ensure fidelity in the process of docking and/or fusion of vesicles with their correct acceptor compartment.

Overview

| | |
|----------------------|--|
| Product Name | Anti-RAB5A/B/C Antibody Picoband® Biotin Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N. |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. |
| Host | Rabbit |
| Uniprot ID | P20339/P61020/P51148 |

Technical Details

| | |
|---------------------|--|
| Immunogen | E.coli-derived human RAB5A/B/C recombinant protein (Position: M1-N215). |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Biotin |
| Suggested Dilutions | The intended application should be selected according to the customer's experimental requirements. |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-RAB5A/B/C Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.