

## Anti-PAR-3/PARD3 Antibody Picoband® APC Conjugated

Catalog Number: A01897-2-APC

### About PARD3

Partitioning defective 3 homolog is a protein that in humans is encoded by the PARD3 gene. This gene encodes a member of the PARD protein family. PARD family members interact with other PARD family members and other proteins; they affect asymmetrical cell division and direct polarized cell growth. Multiple alternatively spliced transcript variants have been described for this gene.

### Overview

|                      |   |
|----------------------|---|
| Product Name         | Anti-PAR-3/PARD3 Antibody Picoband® APC Conjugated  |
| Reactive Species     | Human   |
| Application          | Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs. |
| Clonality            | Polyclonal  |
| Formulation          | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .  |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.  |
| Host                 | Rabbit  |
| Uniprot ID           | Q8TEW0  |

### Technical Details

|                     |  |
|---------------------|--|
| Immunogen           | E.coli-derived human PAR-3/PARD3 recombinant protein (Position: M1-Q1346). |
| Cross Reactivity    | No cross-reactivity with other proteins.                                   |
| Isotype             | Rabbit IgG   |
| Form                | Liquid   |
| Concentration       | 0.5 mg/mL  |
| Purification        | Immunogen affinity purified.   |
| Conjugate           | APC<br>Excitation Wavelength: 633-647 nm<br>Emission Wavelength: 660 nm    |
| Suggested Dilutions | Optimal dilutions should be determined by end users.                       |

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PAR-3/PARD3 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.