

Anti-C5a-R/C5ar1 Antibody Picoband® HRP Conjugated

Catalog Number: A01898-2-HRP

About C5ar1

The C5a receptor also known as complement component 5a receptor 1 (C5AR1) or CD88 (Cluster of Differentiation 88) is a G protein-coupled receptor for C5a. It functions as a complement receptor. C5a receptor modulates inflammatory responses, obesity, development and cancers.

Overview

Product Name	Anti-C5a-R/C5ar1 Antibody Picoband® HRP Conjugated
Reactive Species	Mouse
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P30993

Technical Details

Immunogen	E.coli-derived mouse C5a-R/C5ar1 recombinant protein (Position: H29-D327).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-C5a-R/C5ar1 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.