

Anti-CaMKIV (phospho-T196/200) CAMK4 Antibody

Catalog Number: A01905T196


Overview

Product Name	Anti-CaMKIV (phospho-T196/200) CAMK4 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-CaMKIV (phospho-T196/200) CAMK4 Antibody catalog # A01905T196. Tested in IHC applications. This antibody reacts with Human, Mouse, Rat.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q16566

Technical Details

Immunogen	Synthetic phosphopeptide derived from human CaMKIV around the phosphorylation site of Threonine 200.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200

Anti-CaMKIV (phospho-T196/200) CAMK4 Antibody (A01905T196) Images



Immunohistochemistry (IHC) analyzes of p-CaMKIV (T196/200) pAb in paraffin-embedded human brain tissue.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CaMKIV (phospho-T196/200) CAMK4 Antibody

For Research Use Only. Not for use in diagnostic procedures.