

Anti-WD repeat-containing protein 5 WDR5 Antibody

Catalog Number: A01910-2


Overview

Product Name	Anti-WD repeat-containing protein 5 WDR5 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-WD repeat-containing protein 5 WDR5 Antibody catalog # A01910-2. Tested in WB,IHC,Flow Cytometry applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P61964


Technical Details

Immunogen	Recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:1,000 IHC: 1:50-1:200 FC: 1:50-1:100

Anti-WD repeat-containing protein 5 WDR5 Antibody (A01910-2) Images



Western blot analysis of WDR5 on different cell lysates using anti-WDR5 antibody at 1/1,000 dilution. Positive control:
Lane 1: Hela Lane 2: A549 Lane 3: NIH-3T3 Lane 4: PC-12



Immunohistochemical analysis of paraffin-embedded mouse brain tissue using anti-WDR5 antibody. Counter stained with hematoxylin.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-WD repeat-containing protein 5 WDR5 Antibody

For Research Use Only. Not for use in diagnostic procedures.