

## Anti-Histone deacetylase 7 HDAC7 Antibody

Catalog Number: A01913-2


### Overview

Product Name	Anti-Histone deacetylase 7 HDAC7 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Histone deacetylase 7 HDAC7 Antibody catalog # A01913-2. Tested in WB, Flow Cytometry applications. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8WUI4


### Technical Details

Immunogen	recombinant protein
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:2,000 FC: 1:50-1:100

## Anti-Histone deacetylase 7 HDAC7 Antibody (A01913-2) Images



Western blot analysis of HDAC7 on different lysates using anti-HDAC7 antibody at 1/1,000 dilution. Positive control:  
Lane 1: A549 Lane 2: Human brain



Flow cytometric analysis of K562 cells with HDAC7 antibody at 1/50 dilution (red) compared with an unlabelled control (cells without incubation with primary antibody; black). Alexa Fluor 488-conjugated goat anti rabbit IgG was used as the secondary antibody

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Histone deacetylase 7 HDAC7 Antibody

For Research Use Only. Not for use in diagnostic procedures.