

## Anti-TH (phospho-S8) Antibody

Catalog Number: A01917S8


### Overview

Product Name	Anti-TH (phospho-S8) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TH (phospho-S8) Antibody catalog # A01917S8. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P07101

### Technical Details

Immunogen	Synthetic peptide derived from human TH around the phosphorylation site of S8.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

## Anti-TH (phospho-S8) Antibody (A01917S8) Images



Western blot (WB) analysis of p-TH (S8) pAb at 1:500 dilution Lane1:MCF-7 whole cell lysate(40ug) Lane2:MCF-7 treated with Forskolin(4μg/ml,30 minutes) whole cell lysate(40ug) Lane3:MCF-7 treated with Forskolin(4μg/ml,60 minutes) whole cell lysate(40ug)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TH (phospho-S8) Antibody

For Research Use Only. Not for use in diagnostic procedures.