

Anti-Radixin RDX Antibody

Catalog Number: A01926


Overview

Product Name	Anti-Radixin RDX Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Radixin RDX Antibody catalog # A01926. Tested in WB, ICC/IF, IHC, IP, Flow Cytometry applications. This antibody reacts with Human, Mouse, Rat.
Application	Flow Cytometry, IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P35241


Technical Details

Immunogen	recombinant protein
Isotype	Antiserum
Form	Liquid
Concentration	85 mg/mL by Refractometry
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:1,000-1:2,000 ICC: 1:50-1:200 IHC: 1:50-1:200 FC: 1:50-1:100

Anti-Radixin RDX Antibody (A01926) Images



Western blot analysis of Radixin on different lysates using anti-Radixin antibody at 1/1,000 dilution. Positive control:
Lane 1: Hela Lane 2: A431 Lane 3: HepG2



ICC staining Radixin in HeLa cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Radixin RDX Antibody

For Research Use Only. Not for use in diagnostic procedures.