

Anti-KGF/FGF7 Antibody Picoband®

Catalog Number: A01931-1

About FGF7

Keratinocyte growth factor is a protein that in humans is encoded by the FGF7 gene. The protein encoded by this gene is a member of the fibroblast growth factor (FGF) family. FGF family members possess broad mitogenic and cell survival activities, and are involved in a variety of biological processes, including embryonic development, cell growth, morphogenesis, tissue repair, tumor growth and invasion. This protein is a potent epithelial cell-specific growth factor, whose mitogenic activity is predominantly exhibited in keratinocytes but not in fibroblasts and endothelial cells. Studies of mouse and rat homologs of this gene implicated roles in morphogenesis of epithelium, reepithelialization of wounds, hair development and early lung organogenesis.

Overview


Product Name	Anti-KGF/FGF7 Antibody Picoband®
Reactive Species	Human
Description	Boster Bio Anti-KGF/FGF7 Antibody Picoband® catalog # A01931-1. Tested in ELISA, WB applications. This antibody reacts with Human. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Each vial contains 4mg Trehalose, 0.9mg NaCl and 0.2mg Na ₂ HPO ₄ .
Storage Instructions	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P21781

Technical Details


Immunogen	E.coli-derived human KGF/FGF7 recombinant protein (Position: R54-K178).
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Lyophilized

Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.25-0.5ug/ml, Human ELISA, 0.1-0.5ug/ml, -

Anti-KGF/FGF7 Antibody Picoband® (A01931-1) Images



Western blot analysis of KGF/FGF7 using anti-KGF/FGF7 antibody (A01931-1). Electrophoresis was performed on a 5-20% SDS-PAGE gel at 70V (Stacking gel) / 90V (Resolving gel) for 2-3 hours. The sample well of each lane was loaded with 50ug of sample under reducing conditions. Lane 1: human U20S whole cell lysates. After Electrophoresis, proteins were transferred to a Nitrocellulose membrane at 150mA for 50-90 minutes. Blocked the membrane with 5% Non-fat Milk/ TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-KGF/FGF7 antigen affinity purified polyclonal antibody (Catalog # A01931-1) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an Enhanced Chemiluminescent detection (ECL) kit (Catalog # EK1002) with Tanon 5200 system. A specific band was detected for KGF/FGF7 at approximately 22KD. The expected band size for KGF/FGF7 is at 22KD.



The inductive ability of DPCs. (A) The intact DPC markers ALP, BMP2, and FGF7 were detected by Western blot. (B-F) HaCaT cells were cultured with conditioned medium from DPCs treated with virus expressing Tcf4, Twist1, and siTwist1. Differentiation markers KRT40, MSX2, KRT5, and KRT15 were detected by Western blot. (C-F) Quantitative analysis of (B). The expression levels of target genes were standardized to that of GAPDH. * P < 0.05 when compared with the control group. ** P < 0.01 when compared with the control group. *** P < 0.001 when compared with the control group. Index in PubMed under a CC BY license. PMID: 32974352

1 Publications Citing This Product

1. PubMed ID: 30823849, Li X,Wang X,Bai L,Zhao P,Zhang M.Exposure to 50 Hz electromagnetic fields enhances hair follicle regrowth in C57BL/6 mice.Exp Biol Med(Maywood).2019 Apr;244(5):389-394.doi:10.1177/1535370219834639.Epub 2019 Mar 1.PMID:30823849;PMCID:PMC6488868.

Visit bosterbio.com/anti-kgf-fgf7-picoband-trade-antibody-a01931-1-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



For Research Use Only. Not for use in diagnostic procedures.