

Anti-HtrA2/Omi Antibody

Catalog Number: A01941-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-HtrA2/Omi Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-HtrA2/Omi Antibody (Catalog# A01941-1). Tested in IHC, ICC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, ICC
Clonality	Polyclonal
Formulation	PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short-term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O43464

Technical Details

Immunogen	A peptide derived from human HtrA2. Immunogen sequence location: 181-230
Cross Reactivity	No cross reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.

Suggested Dilutions


ICC/IF: 1:100 - 1:300.

ELISA: 1:40000.

Not yet tested in other applications.

Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users.

Anti-HtrA2/Omi Antibody (A01941-1) Images



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using HtrA2 Antibody. The picture on the right is blocked with the synthesized peptide.

Anti-HtrA2/Omi Antibody

For Research Use Only. Not for use in diagnostic procedures.