

Anti-Actinin Alpha 4 ACTN4 Antibody

Catalog Number: A01975


Overview

Product Name	Anti-Actinin Alpha 4 ACTN4 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Actinin Alpha 4 ACTN4 Antibody catalog # A01975. Tested in WB,ICC,IF,IP,IHC,Flow Cytometry applications. This antibody reacts with Human,Mouse,Rat.
Application	Flow Cytometry, IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O43707


Technical Details

Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 ICC: 1:100-1:500 IHC: 1:50-1:200 FC: 1:50-1:100 IP: 1:10-1:50

Anti-Actinin Alpha 4 ACTN4 Antibody (A01975) Images



Western blot analysis of alpha Actinin 4 on different lysates using anti-alpha Actinin 4 antibody at 1/500 dilution. Positive control: Lane 1: Hela Lane 2: PC-12 Lane 3: NIH-3T3 Lane 4: Rat liver tissue Lane 5: A431 Lane 6: HepG2



ICC staining alpha Actinin 4 in A549 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Actinin Alpha 4 ACTN4 Antibody

For Research Use Only. Not for use in diagnostic procedures.