

Anti-RRM2 Antibody

Catalog Number: A01978


Overview

Product Name	Anti-RRM2 Antibody
Reactive Species	Human
Description	Boster Bio Anti-RRM2 Antibody catalog # A01978. Tested in WB,IP,ICC/IF applications. This antibody reacts with Human.
Application	IP, IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P31350

Technical Details

Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 ICC: 1:20-1:50 IP: 1:10-1:50

Anti-RRM2 Antibody (A01978) Images



Western blot analysis of RRM2 on HeLa (1) and HepG2 (2) cell lysate using anti-RRM2 antibody at 1/1,000 dilution.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RRM2 Antibody

For Research Use Only. Not for use in diagnostic procedures.