

Anti-AChR alpha3 CHRNA3 Antibody

Catalog Number: A01981


Overview

Product Name	Anti-AChR alpha3 CHRNA3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-AChR alpha3 CHRNA3 Antibody catalog # A01981. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P32297


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human AChRalpha3. AA range:90-139
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

Anti-AChR alpha3 CHRNA3 Antibody (A01981) Images



Western Blot analysis of various cells using AChRalpha3 Polyclonal Antibody



Western blot analysis of lysate from HeLa cells, using AChRalpha3 antibody.

1 Publications Citing This Product

1. PubMed ID: 10.1016/j.ejphar.2019.172674, Identification of nicotinic acetylcholine receptor subunits in different lung cancer cell lines and the inhibitory effect of alpha-conotoxin Tx1D on lung cancer cell growth

Visit bosterbio.com/anti-achr-alpha3-antibody-a01981-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AChR alpha3 CHRNA3 Antibody

For Research Use Only. Not for use in diagnostic procedures.