

Anti-Human CHRNA3 DyLight® 488 conjugated Antibody

Catalog Number: A01981-Dyl488

Overview

Product Name	Anti-Human CHRNA3 DyLight® 488 conjugated Antibody
Reactive Species	Human
Description	Boster Bio Anti-Human CHRNA3 DyLight® 488 conjugated Antibody catalog # A01981-Dyl488. Tested in Flow Cytometry applications. This antibody reacts with Human.
Conjugate	DyLight®488
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P32297

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human CHRNA3, which shares 97.4% and 94.7% amino acid (aa) sequence identity with mouse and rat CHRNA3, respectively.
Predicted Reactive Species	Human
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used: Flow Cytometry, 1-3ug/1x10 ⁶ cells

Anti-Human CHRNA3 DyLight® 488 conjugated Antibody (A01981-Dyl488) Images



Boster Kit Box

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Human CHRNA3 DyLight® 488 conjugated Antibody