

## Anti-Phospho-Gab 1 (Y307) Antibody

Catalog Number: A01989Y307

### Overview


Product Name	Anti-Phospho-Gab 1 (Y307) Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Phospho-Gab 1 (Y307) Antibody catalog # A01989Y307. Tested in IHC, IF, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, IHC
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q13480

### Technical Details

Immunogen	Synthesized phospho-peptide around the phosphorylation site of human Gab 1 (phospho Tyr307)
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	IHC 1:100-1:300 ELISA 1:40000 IF 1:50-200

## Anti-Phospho-Gab 1 (Y307) Antibody (A01989Y307) Images

---



Immunohistochemical analysis of paraffin-embedded human Colon cancer. 1, Antibody was diluted at 1:200 (4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200 (room temperature, 45min).

### Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Phospho-Gab 1 (Y307) Antibody

For Research Use Only. Not for use in diagnostic procedures.