

Anti-Lsk MATK Monoclonal Antibody

Catalog Number: A02050


Overview

Product Name	Anti-Lsk MATK Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-Lsk MATK Monoclonal Antibody catalog # A02050. Tested in WB, Flow Cytometry, ELISA applications. This antibody reacts with Human.
Application	ELISA, Flow Cytometry, WB
Clonality	Monoclonal 9D7
Formulation	Antibody are purified by protein G affinity chromatography. Liquid in PBS containing 50% glycerol and 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P42679


Technical Details

Immunogen	Purified recombinant fragment of human Lsk expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Protein A purified
Suggested Dilutions	WB 1:500-1:2000 FC 1:200-1:400 ELISA 1:10000

Anti-Lsk MATK Monoclonal Antibody (A02050) Images



Western Blot analysis using Lsk Monoclonal Antibody against K562 cell lysate (1).



Flow cytometric analysis of K562 cells using Lsk Monoclonal Antibody (green) and negative control (purple).

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Lsk MATK Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.