

Anti-Fnk PLK3 Antibody

Catalog Number: A02054-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-Fnk PLK3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Fnk PLK3 Antibody catalog # A02054-1. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9H4B4


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human PLK3. AA range:231-280
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000


ELISA 1:10000

This antiserum has been tested for use in western blotting, immunoprecipitation and immunohistochemistry. Specific conditions for reactivity should be optimized by the end user. Expect bands approximately 150kDa by western blotting in the appropriate cell lysate or extract.


Anti-Fnk PLK3 Antibody (A02054-1) Images



Western blot analysis of the lysates from RAW264.7 cells using PLK3 antibody.



Western Blot analysis of various cells using Fnk Polyclonal Antibody



Western blot analysis of lysates from HepG2 cells, using PLK3 Antibody. The lane on the right is blocked with the synthesized peptide.

Anti-Fnk PLK3 Antibody

For Research Use Only. Not for use in diagnostic procedures.