

Anti-Glutathione Peroxidase 4/GPX4 Antibody Picoband®

Catalog Number: A02059-1

About GPX4

Glutathione peroxidase 4, also known as GPX4, is an enzyme that in humans is encoded by the GPX4 gene. This gene encodes a member of the glutathione peroxidase protein family. Glutathione peroxidase catalyzes the reduction of hydrogen peroxide, organic hydroperoxide, and lipid peroxides by reduced glutathione and functions in the protection of cells against oxidative damage. Human plasma glutathione peroxidase has been shown to be a selenium-containing enzyme and the UGA codon is translated into a selenocysteine. The encoded protein has been identified as a moonlighting protein based on its ability to serve dual functions as a peroxidase as well as a structural protein in mature spermatozoa. Through alternative splicing and transcription initiation, rat produces proteins that localize to the nucleus, mitochondrion, and cytoplasm. In humans, alternative transcription initiation and the cleavage sites of the mitochondrial and nuclear transit peptides need to be experimentally verified. Alternative splicing results in multiple transcript variants.


Overview

Product Name	Anti-Glutathione Peroxidase 4/GPX4 Antibody Picoband®
Reactive Species	Human, Monkey, Mouse, Rat
Description	Boster Bio Anti-Glutathione Peroxidase 4/GPX4 Antibody Picoband® catalog # A02059-1. Tested in ELISA, Flow Cytometry, IHC, WB applications. This antibody reacts with Human, Monkey, Mouse, Rat. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Each vial contains 4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P36969


Technical Details

Immunogen	E.coli-derived human Glutathione Peroxidase 4/GPX4 recombinant protein (Position: A30-F197).
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG


Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.25-0.5 ug/ml, Human, Monkey, Mouse, Rat Immunohistochemistry(Paraffin-embedded Section), 2-5 ug/ml, Human, Mouse Flow Cytometry(Fixed), 1-3 ug/1x10 ⁶ cells, Human ELISA, 0.1-0.5 ug/ml, -




Western blot analysis of Glutathione Peroxidase 4/GPX4 using anti-Glutathione Peroxidase 4/GPX4 antibody. Electrophoresis was performed on a 12% SDS-PAGE gel at 80V (Stacking gel) / 120V (Resolving gel) for 4 h. One well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: human HepG2 whole cell lysates, Lane 2: human MCF-7 whole cell lysates, Lane 3: human RT4 whole cell lysates, Lane 4: monkey COS7 whole cell lysates, Lane 5: rat kidney tissue lysates, Lane 6: rat testis tissue lysates, Lane 7: mouse kidney tissue lysates. The proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-Glutathione Peroxidase 4/GPX4 antibody affinity purified polyclonal antibody (Catalog # A02059-1) at 0.5 ug/mL overnight at 4°C, then washed with TBST 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:1000 for 1 hour at RT. The signal is developed using an ECL Plus Western Blotting Substrate (Catalog # AR1196) system. A specific band was detected for Glutathione Peroxidase 4/GPX4 at approximately 19 kDa. The molecular weight for Glutathione Peroxidase 4/GPX4 is at 22 kDa.



Effects of MCTR1 on ferroptosis in CLP-induced AKI. Mice were given MCTR1 once, twice, or not during CLP. The once-daily administration mode (qd): Mice were given MCTR1 (200 ng/mice, iv) 0.5 h before being subjected to CLP. The twice-daily administration mode (bid): Mice were given MCTR1 (200 ng/mice, iv) 0.5 h before being subjected to CLP, and an additional MCTR1 (200 ng/mice, iv) was given 12 h after CLP. All kidney samples were collected at 24 h after CLP. Shown are representative western blotting and quantification of GPX4 and PTGS2. C-E Quantitative analyses of MDA, GSH, and non-heme iron. F Quantitative analyses of Perl's stain. G Representative Perl's stain images. H Representative TEM images. The black arrow indicates ferroptosis-like mitochondria. I Quantitative analyses of mitochondrial length. J Quantitative analyses of renal tubular injury. K, L Shown are representative HE stains and PAS stain. M Representative IHC images. N Quantitative analyses of IHC stain of NGAL. O, P Quantitative analyses of serum Scr and BUN. n = 6 mice/group. **P < 0.01. iv: intravenous Index in PubMed under a CC BY license. PMID: 34961563



Effects of Nrf2 on MCTR1-regulated ferroptosis in CLP-induced AKI. Mice were given ML385 (30 mg/kg, iv) 0.5 h before being subjected to CLP, and a twice-daily administration mode of MCTR1 as described before in the text. All kidney samples were collected at 24 h after CLP. A-D Shown are representative western blotting and quantification of GPX4 and PTGS2. E - G Quantitative analyses of MDA, GSH, and non-heme iron. H Quantitative analyses of Perl's stain. I Representative Perl's stain images. J Representative TEM images. The black arrow indicates ferroptosis-like mitochondria. K Quantitative analyses of mitochondrial length. L Histological analyses of renal tubular injury. M, N Representative HE and PAS stain. O Representative IHC images for NGAL. P Quantitative analyses of IHC stain of NGAL. Q Quantitative analyses of serum Scr and BUN. n = 6 mice/group, mean ± SD were presented. **P < 0.01 Index in PubMed under a CC BY license. PMID: 34961563



Effects of MCTR1 on LPS-induced lipid peroxidation and ferroptosis in vitro. A Viability curves for HK-2 cells treated with different concentrations (0, 1, 3, 10, 30, 100, and 300 nM) of MCTR1 and LPS (1 ug/ml) for 8 h. B Viability curves for HK-2 cells treated with different concentrations (0, 1, 3, 10, 30, 100, and 300 nM) of MCTR1 for 0 and 8 h. C HK-2 cells treated with MCTR1 (100 nM) and LPS (1 ug/ml) for 8 h, and the visualization of cell viability were evaluated by phase-contrast microscopy. D, E Representative Western blotting and quantification of GPX4 and PTGS2 protein. F - H Quantitative analyses of MDA, GSH, and non-heme iron. I Quantitative analyses of oxidized C11-BODIPY 581/591 probe by flow cytometry. J Representative flow cytometry images of C11-BODIPY 581/591 fluorescent ratio-probe. K Representative TEM images. The black arrow indicates ferroptosis-like mitochondria. L Quantitative analyses of mitochondrial length. n = 3, mean ± SD were presented. Index in PubMed under a CC BY license. PMID: 34961563

peroxidase (PHGPx)

2. PubMed ID: 10.1186/s13578-021-00734-x, Maresin conjugates in tissue regeneration-1 suppresses ferroptosis in se

Visit bosterbio.com/anti-glutathione-peroxidase-4-gpx4-picoband-trade-antibody-a02059-1-boster.html to see all 2 p

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Glutathione Peroxidase 4/GPX4 Antibody

For Research Use Only. Not for use in diagnostic procedures.