

Anti-Deleted in azoospermia-like DAZL Antibody

Catalog Number: A02069


Overview

Product Name	Anti-Deleted in azoospermia-like DAZL Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Deleted in azoospermia-like DAZL Antibody catalog # A02069. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q92904


Technical Details

Immunogen	recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500 IHC: 1:50

Anti-Deleted in azoospermia-like DAZL Antibody (A02069) Images



Immunohistochemical analysis of paraffin-embedded rat testes tissue using anti-DAZL antibody. Counter stained with hematoxylin.



Immunohistochemical analysis of paraffin-embedded mouse testes tissue using anti-DAZL antibody. Counter stained with hematoxylin.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Deleted in azoospermia-like DAZL Antibody

For Research Use Only. Not for use in diagnostic procedures.