

Anti-Fumarylacetoacetase FAH Antibody

Catalog Number: A02072


Overview

Product Name	Anti-Fumarylacetoacetase FAH Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Fumarylacetoacetase FAH Antibody catalog # A02072. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P16930


Technical Details

Immunogen	Recombinant full length Human FAH.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:50-1:200

Anti-Fumarylacetoacetase FAH Antibody (A02072) Images



WesternBlot (WB) analysis of FAH polyclonal antibody



WesternBlot (WB) analysis of FAH polyclonal antibody

1 Publications Citing This Product

1. PubMed ID: 10.1016/j.biopha.2018.11.148, Urinary proteome analysis of acute hypercoagulable state in rat model induced by epsilon-aminocaproic acid

Visit bosterbio.com/anti-fumarylacetoacetase-fah-antibody-a02072-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Fumarylacetoacetase FAH Antibody

For Research Use Only. Not for use in diagnostic procedures.