

Anti-TRIM21 Antibody

Catalog Number: A02079-1


Overview

Product Name	Anti-TRIM21 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TRIM21 Antibody catalog # A02079-1. Tested in WB,IHC,ICC/IF,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P19474


Technical Details

Immunogen	Recombinant fusion protein of human TRIM21/SS-A(NP_003132.2).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:50-1:200 ICC/IF: 1:50-1:200 IP: 1:50-1:200


Anti-TRIM21 Antibody (A02079-1) Images



Western blot analysis of extracts of various cell lines, using TRIM21/SS-A antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 3s.



Immunohistochemistry of paraffin-embedded mouse spinal cord using TRIM21/SS-A Rabbit pAb at dilution of 1:50. Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.



Immunohistochemistry of paraffin-embedded rat kidney using TRIM21/SS-A Rabbit pAb at dilution of 1:50. Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.

1 Publications Citing This Product

1. PubMed ID: -, Chun Wang, Bin Wang, Maimaiti Aili, Shixia Huo, Meng-ting Han, Silafu Aibai, Zhi-Jian Li, "Effect of Artemisia rupestris L. Extract on Gastrointestinal Hormones and Brain-Gut Peptides in Functional Dyspepsia Rats", Evidence-Based Complementary and Alternative M

Visit bosterbio.com/anti-trim21-antibody-a02079-1-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TRIM21 Antibody

For Research Use Only. Not for use in diagnostic procedures.