

## Anti-Prostatic acid phosphatase ACPD Antibody

Catalog Number: A02082-1


### Overview

Product Name	Anti-Prostatic acid phosphatase ACPD Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Prostatic acid phosphatase ACPD Antibody catalog # A02082-1. Tested in WB,IHC,ICC/IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P15309


### Technical Details

Immunogen	Recombinant fusion protein of human Prostatic Acid Phosphatase (PAP/Prostatic Acid Phosphatase (ACPD))(NP_001090.2).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:50-1:200 ICC/IF: 1:50-1:200


## Anti-Prostatic acid phosphatase ACP (A02082-1) Images




Western blot analysis of extracts of various cell lines, using Prostatic Acid Phosphatase antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 180s.



Immunohistochemistry of paraffin-embedded Human benign prostatic hyperplasia using Prostatic Acid Phosphatase Rabbit pAb at dilution of 1:50. Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.



Immunohistochemistry of paraffin-embedded Human colon using Prostatic Acid Phosphatase Rabbit pAb at dilution of 1:50. Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.



Immunohistochemistry of paraffin-embedded human kidney using Prostatic Acid Phosphatase Rabbit pAb at dilution of 1:50. Perform high pressure antigen retrieval with 10 mM citrate buffer pH 6.0 before commencing with IHC staining protocol.

**Submit a product review to Biocompare.com**

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Prostatic acid phosphatase ACPA Antibody

For Research Use Only. Not for use in diagnostic procedures.