

Anti-Exostosin-1 EXT1 Antibody

Catalog Number: A02086

Overview

Product Name	Anti-Exostosin-1 EXT1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Exostosin-1 EXT1 Antibody catalog # A02086. Tested in IHC,IF,ICC,ELISA applications. This antibody reacts with Human,Mouse,Rat.
Application	ELISA, IF, IHC, ICC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q16394

Technical Details

Immunogen	KLH conjugated synthetic peptide derived from human EXT1:651-746/746
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC : 1:100-1:500 IF: 1:100-1:500 ICC: 1:100-1:500 ELISA: 1:5000-1:10000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Exostosin-1 EXT1 Antibody

For Research Use Only. Not for use in diagnostic procedures.