

## Anti-HIF1AN Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A02116-2-Fluoro488

### About HIF1AN

Enables several functions, including 2-oxoglutarate-dependent dioxygenase activity; NF-kappaB binding activity; and transition metal ion binding activity. Involved in negative regulation of Notch signaling pathway and positive regulation of myoblast differentiation. Located in cytosol; nucleoplasm; and perinuclear region of cytoplasm.

### Overview

Product Name	Anti-HIF1AN Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9NWT6

### Technical Details

Immunogen	E.coli-derived human HIF1AN recombinant protein (Position: D28-N349). Human HIF1AN shares 98.4% amino acid (aa) sequence identity with mouse HIF1AN.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-HIF1AN Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.