

Anti-LGR4 Antibody

Catalog Number: A02134-1

About LGR4

LGR4, also known as GPR48, is a membrane of the leucine-rich GPCR family. LGR4 is composed of a large NH₂-terminal extracellular domain containing 18 LRRs, the seven membrane-spanning segments, and a COOH-terminal intracellular domain. LGR4 is a receptor for R-spondins that potentiates the canonical Wnt signaling pathway and is involved in the formation of various organs. It has also been shown that LGR4 regulates cyclic AMP/protein kinase A (PKA) signalling. (PMID: 23756652; 17178856; 24212090)

Overview

Product Name	Anti-LGR4 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-LGR4 Antibody catalog # A02134-1. Tested in IHC, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IHC
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9BXB1

Technical Details

Immunogen	E.coli-derived human LGR4 recombinant protein (Position: C29-S373).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Immunohistochemistry, 1:50-400 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LGR4 Antibody

For Research Use Only. Not for use in diagnostic procedures.