

## Anti-Nuclear RNA export factor 1 NXF1 Antibody

Catalog Number: A02137-1


### Overview

Product Name	Anti-Nuclear RNA export factor 1 NXF1 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Nuclear RNA export factor 1 NXF1 Antibody catalog # A02137-1. Tested in WB,ICC,IP applications. This antibody reacts with Human.
Application	IP, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UBU9


### Technical Details

Immunogen	Recombinant protein within c-terminal human NXF1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 ICC: 1:50-1:100 IP: 1:10-1:50

## Anti-Nuclear RNA export factor 1 NXF1 Antibody (A02137-1) Images



Western blot analysis of NXF1 on 293 (1) and SiHa (2) cell lysate using anti-NXF1 antibody at 1/2,000 dilution.



ICC staining NXF1 in LOVO cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Nuclear RNA export factor 1 NXF1 Antibody

For Research Use Only. Not for use in diagnostic procedures.