

Anti-AOX1/Aldehyde Oxidase Antibody

Catalog Number: A02144-1


Overview

Product Name	Anti-AOX1/Aldehyde Oxidase Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-AOX1/Aldehyde Oxidase Antibody catalog # A02144-1. Tested in WB applications. This antibody reacts with Human, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q06278

Technical Details

Immunogen	Synthesized peptide derived from the Internal region of human AOX1. at AA range: 490-570
Predicted Reactive Species	Bovine, Canine, Equine, Pig, Rabbit
Isotype	IgG
Form	Liquid
Concentration	This antibody's concentration is >0.5mg/ml.
Purification	Immunogen affinity purified
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used: WB, 1:500-1:2000

Anti-AOX1/Aldehyde Oxidase Antibody (A02144-1) Images



Western Blot (WB) analysis of specific cells using AOX1 Polyclonal antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-AOX1/Aldehyde Oxidase Antibody